双极电极
Bipolar Electrodes

双极电极
Bipolar Electrodes

REF:
700450 – 700499

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STERRAD® is a trademark of Johnson & Johnson, Inc.

English

EN

Product / User / Disposal:
Electrosurgery accessories may only be used and disposed of by qualified medical personnel!
Notwithstanding these instructions, reading the instructions for use provided with the electrosurgery device and other accessories being used is mandatory.
△ **Not sterile.** Must be cleaned and sterilized prior to initial and each subsequent use.

Intended use:
Bipolar electrodes for the coagulation / volume reduction / cutting of soft tissue during surgical procedures.

Service life:
A minimum of 20 reconditioning cycles can be expected with proper use.

Before using:
△ Inspect the product for cleanliness, intact insulation, and integrity before each use. Note the sensitive instrument tips in particular.
△ Only use sterilized products that are in flawless condition! All products must be used with a suitable connection cable from Sutter!
△ To connect electrodes or accessories, the electrosurgery device must be turned off or in standby mode. Failure to comply may lead to burns and electric shocks!

During use:
Operating guidelines for these product groups are available as non-binding instructions. They can be requested from the manufacturer.
△ Always work with the lowest output setting for the desired surgical effect.
△ Electrode tips may cause injuries!
△ After application, electrode tips may be so hot they cause burns!
△ Never lay down the electrode on or in the immediate vicinity of the patient! Lay cables and store unused instruments so they are isolated from the patient.
△ Not for use in the presence of combustible or explosive substances!
△ Unintentional activation or movement of the activated electrode outside the field of vision can lead to unwanted burns!
△ Maximum allowable voltage 500 Vp.

Reconditioning:
Observe national guidelines and regulations!
Disconnect electrodes from accessories / connection cables!
The overall reconditioning process encompasses pre-cleaning, cleaning / disinfection, and sterilization.
△ Machine cleaning / disinfection is preferred due to effectiveness and reproducibility!
△ Do not immerse in hydrogen peroxide (H₂O₂)!

Pre-cleaning:
Do not allow blood and tissue residues to dry up; rinse thoroughly with cold water after max. 1 hour! Use a soft brush if needed (no wire brush or similar).

Manual cleaning and disinfection:	
Cleaning step	Description
Pre-cleaning	Rinse for 5 minutes under cold water and clean with a soft brush (such as MED100.33 from Medisafe GmbH) until no more residues are visible.
Ultrasound and disinfection	Ultrasound bath 35 kHz at room temperature, 10 minutes, cleaning/ disinfection solution 2 % Bomix® plus (Bode Chemie).
Secondary cleaning	Rinse hard to reach areas for 20 seconds with a spray nozzle as needed, then rinse the entire instrument for 30 seconds with demineralized water.

Proof of fundamental suitability of the instruments for manual cleaning and disinfection was provided for a worst-case product by an independent, accredited test laboratory using the instrument disinfectant Bomix® plus (Bode Chemie) with a concentration of 2 %, test report no. 07015-2 dated 24.11.2015. Transferability was proven by the internal validation no. V278.

Machine cleaning and disinfection:
In selecting the cleaning and disinfection device (CDD), note that the effectiveness of the CDD must have been tested (DGHM ("Deutsche Gesellschaft für Hygiene und Mikrobiologie e.V.", German Association for Hygiene and Microbiology) or FDA clearance for instance, or the CE marking according to EN ISO 15883).
• Load instruments into the CDD. In doing so, make sure the instruments do not touch each other and are securely supported.

Program steps	Parameters
Pre-rinse	10±2 °C, 1 minute
Cleaning with 0.3 % (3 ml/liter) deconex® 28 ALKA ONE-x	70±2 °C, 5 minutes
Final rinse	10±2 °C, 2 minutes
Thermal disinfection	90±2 °C, 5 minutes

Proof of fundamental suitability of the instruments for effective machine cleaning and disinfection was provided by an independent, accredited test laboratory for a worst-case product using the Miele G7836 CD disinfecter (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh, Germany) and the cleaning agent deconex® 28 ALKA ONE-x with a concentration of 0.3 % (Borer Chemie AG, Zuchwil, Switzerland), Test report no. 100331-10 dated 30.04.2010. Transferability was proven by the internal validation no. V283.

• Please note: The preceding instructions contain validated minimum times for successful cleaning with the prescribed program steps. Deviating process parameters (longer cleaning time and higher cleaning temperatures up to 95 °C) do not damage the instruments and are permitted according to the A₀ concept, comparable A₀-value>3000. If a different cleaning agent is used, only use products with properties comparable to deconex® 28 ALKA ONE-x (Borer Chemie), for instance regarding the pH value and compatibility with plastics. Please contact your supplier or hygiene officer in case of questions.

Inspection:
Prior to sterilization, visually inspect the instrument and check for intact insulation, cleanliness, and integrity.

Maintenance:
None

Packaging:
Package cleaned and disinfected instruments in disposable sterilization packaging (single or double packaging – note sharp instrument tips!), or store the instrument in a suitable sterilization container meeting the following requirements:
• EN ISO/ANSI AAMI ISO 11607
• Suitable for steam sterilization (resistant to temperatures up to min. 141 °C, adequate vapor permeability)
• Sufficient protection of the instruments and/or sterilization packaging against mechanical damage

Sterilization:
Products must be cleaned and disinfected prior to sterilization.
• Steam sterilization, steam sterilizer according to EN 13060 and/or EN 285 and validated according to EN ISO 17665

Program steps	Parameters
Process	Fractionated vacuum (dynamic evacuation)
Sterilization temperature	132 °C (max. 138 °C plus tolerance according to EN ISO 17665)
Sterilization time (holding time at sterilization temperature)	min. 3 minutes
Drying time	30 minutes

Proof of fundamental suitability of the instruments for effective steam sterilization was provided for a worst-case product by an independent, accredited test laboratory, test report no. 100332-10 dated 19.05.2010. Typical conditions in clinics and medical practices as well as the process described above were taken into account. Transferability was proven by the internal validation no. V283.

△ Do not sterilize in hot air!
△ Do not sterilize in STERRAD®!
△ In case of potential contact with prions (CJD – risk of contamination), destroy the instrument and do not use it again.

Storage / Transportation:
Store in a dry place. Protect from sunlight. Store and transport in secure containers / packaging.
For returns, products must be cleaned, disinfected, and packed in sterile packaging.

Please note:
Any change to the product or deviation from these instructions for use waives the liability of Sutter Medizintechnik GmbH.
Changes reserved.

STERRAD® is a trademark of Johnson & Johnson, Inc.